



Cytology (FNA) Samples



FNA Slides

>5 Unstained + 1 Diffquick



FNA Liquid

>0.5-mL in additive-free tube

Histology (Biopsy) Samples



FFPE Scrolls

10 x 10-micron scrolls
+ 1 H&E Slide



FFPE Slides

10 unstained
+ 1 H&E



Paraffin Block



Fresh Frozen Tissue

6–7 x 30-mg pieces

Non-decalcified



bone samples

Vidium can manage the transport of samples from **YOUR LAB TO OURS**